

Thanks to the use of the latest technologies, the **SECURE** permanent make-up cartridges allow for unprecedented precision. There is definitely more pigment left in the skin and the treatment is practically painless.

S

Ш

(1)

 \cap

 \square

 \square

1

 \Box

Ω

Ш

4

Σ

Z

Ш

Z

4

Σ

 $\mathbf{\alpha}$

Ш

Ω

SECURE cartridges are the first ones in the world that are produced in two versions: with acupuncture needles and traditional needles. Thanks to the innovative use of acupuncture needles in cartridges, using the acupuncture technique is now not only extremely precise, but also safe.

The **PROTECT** technology (patent pending) keeps the needle sterile clean and prevents the dye and blood from returning to the handpiece. Thus, it eliminates the possibility of infection during the treatment. Working with **SECURE** cartridges meets all the requirements of future legal regulations announced by the Ministry of Health regarding the safety of permanent make-up.



Work on the **SECURE** cartridge lasted nearly 5 years. This time was needed for testing the best solutions in the field of: needle blade shape, sharpening technology, needle hardness / elasticity, dye flow, hygiene barrier, spring element mechanics and construction of the production line. The highest quality and precision are our priorities. Thanks to the use of the latest technologies, we have achieved a level that allows the production of individual parts of the cartridge with an accuracy of up to 0.01 millimetres. This translates into the highest quality and precision of work.

The **SECURE** cartridges have been developed entirely by Polish engineers and are manufactured and sterilized in Poland.



CARTRIDGES

ACUPUNCTURE NEEDLES



TRADITIONAL NEEDLES



Innovative cartridges with acupuncture needles and the **PROTECT** system ensure the highest precision as well as hygiene and safety of the treatment.

Extremely sharp, flexible needle gently injects the skin without jerking it or causing excessive irritation.

Bleeding is kept to a minimum, so that much more pigment remains in the skin. As a result, the procedure is practically painless. The skin heals very quickly following the treatment.




ACUPUNCTURE **NEEDLES**

- 1 liner 0.14 HR
- 1 liner 0.16 HR
- 1 liner 0.18 HR
- 1 liner 0.20 HR 1 liner 0.23 HR
- 1 liner 0.25 HR
- 1 liner 0.30 HR
- 1 liner 0.35 HR
- 1 liner 0.40 HR

bioEvolution





90°, 60°, 45°

• short, medium, long needle ejection.

• short ejection of the needle

• stretched out skin, delicate injection.



Cartridges with traditional needles, thanks to the use of stiff needles are ideal for pigmenting both delicate and hard skins.

Thanks to the properly selected hardness of the needle the injection is very precise, so we can minimize skin injuries and bleedina.

The quick and strong entry of the needle deep into the skin means that we get an even, smooth line of injections. As a result of using such solutions, more pigment is left in the skin and the healing time following the treatment is significantly reduced.


TRADITIONAL NEEDLES

1 liner	0.20
• 1 liner	0.25
1 liner	0.30
• 1 liner	0.35
• 1 liner	0.40
3 liner	0.30
3 micro	0.18
3 outline	0.25
• 3 slope	0.40
4 flat	0.35
5 round	0.30
7 round	0.30
• 7 flat convex	0.20

• available soon

Meso





The INFINITIA permanent make-up handpiece has been designed for the most demanding linergists. The innovative drive system ensures the highest precision of the needle. Thanks to the modern and ergonomic design - even long work with the handpiece is not tiring. The highest quality anodized aluminium was used for the production of the device. Thanks to the properly selected weight and excellent workmanship the handpiece operation is very quiet.



Distributor:

www.bioevolution.com